BR1
Fintek Bio-Electric Inc.

Pulsed Electromagnetic Field, Capacitively Coupled Bone Growth Stimulator
Health Canada Class II Medical Device
Manufactured in Canada by:
Fintek-Bio Electric INC.
An ISO13485 company
344 Townsend Avenue
Burlington, Ontario L7T 2A4
1-866-440-5692
info@fintek.ca
www.fintek.ca

Device Description

- The Bio-Electric Response is a Pulsed Electromagnetic Field Capacitive Coupled device that delivers electrical impulses at a range of low frequency known to stimulate the healing process of fractured bones and surrounding damaged tissue. This specific frequency level has been found effective in healing bone fractures that would not normally heal, particularly in patients that may have additional health or injury factors that increase the risk of delayed/non-union.
- The components of the device are housed in a lightweight plastic casing making it a portable hand-held device. The casing clips to your belt, or can be kept in your pocket. The unit runs on a standard 9V DC alkaline battery, which provides over 50 hours of use.
- The device uses reusable, self-adhering electrodes which are generally placed on opposing sides of the fracture site utilizing an alternating current.
- The lead wire that connects the Bio-Electric Response to the electrodes is a 3.5 mm mono plug jack at the device end with two strands and .080 plugs that connect to the electrodes.

Contents of Each Package:

- One Bio-Electric Response Unit
- Two sets of leads
- 9 – 9V Duracell batteries
- 4 – Pairs of Electrodes

Instructions for Proper Use

In order for the Bio-Electric Response to perform its intended use, the follow procedure must be followed:

1. Connect Electrodes to the lead wire

2. Connect the lead wire to the jack located on the Fintek Bio-Electric Response.

3. Remove the plastic pads on the electrodes; place the electrodes directly on either side of the fracture site over the skin to create a circuit.

4. Move the power switch to the on position. The flashing green light indicates a circuit has been formed and treatment is in process. Keep the electrode pads placed on either side of the injury site until full treatment time has elapsed (three hours). The device will stop treatment automatically after the three hours has elapsed, but you must ensure to turn off the unit by switching the power button to "off". This will help to prolong the life of the device's battery which will continue to drain if the switch is left in the "ON" position.
5. Wear the device every day consistently until the fracture has healed or your doctor has notified you that your treatments are complete.
6. Place the electrodes back on the plastic pads in the plastic bag to extend the life of electrodes.

Cautions for Use

The Bio-Electric Response should not be utilized by patients that are pregnant, skeletally immature, or have an implanted cardiac pacemaker without consulting a physician.

Contraindications

The Bio-Electric Response unit should not be used by patients with synovial pseudoarthrosis under any circumstances.

Low Battery and the Replacement of Battery Instructions

There is a small red LED light located on the response unit. When the light turns and stays red and you hear one long beep, this indicates that the battery is running low and should be replaced. Before replacing the used battery, make sure that the unit is turned off. To replace the used battery, open the compartment housing the 9V DC battery, remove the used battery and replace with a new one, observing the correct polarity as shown in the diagram. It is possible to damage the unit by attempting to operate it in reverse polarity. For battery removal, please follow your local guidelines for disposing of batteries.

Proper Storage Conditions

The Bio-Electric Response Bone Growth Stimulator should be kept in a dry place at room temperature. It should not be kept in an area with high humidity as this could cause damage to the device's transformer. The device should never come in contact with water as this could also cause damage to the unit.

Placement of Electrodes

In order for the Response unit to perform its intended use, the electrodes should be placed on opposing sides of the fracture site to create an arc. Ensure that the area where electrodes are to be placed is clean and dry to maintain the life of the electrodes. In cases where the electrodes must be placed on the same surface, i.e. spine, the electrodes need to be placed at least four inches apart to ensure an effective circuit is established with the bone.

Advantages of PEMF

The Fintek Bio-Electric Response employs Pulsed Electromagnetic Fields (PEMF) through capacitively coupled current to promote the healing of bone fractures. The signal targets a bone cell protein, serving as a cell receptor trigger to initiate the proliferation of new bone cells called osteoblasts. The osteoblasts are the building blocks of new bone and effectively fill in any gaps patients may have in their fractures. As an additional benefit, PEMF increases blood flow or vascularization to the treated area, potentially also healing injured tissues, ligaments, tendons, and nerves. The Bio-Electric Response can be used in both appendicular and spinal applications. The device can be prescribed for the following fracture types:

- Open and Closed Fractures
- Complex and Comminuted Fractures
- Spiral, Transverse and Oblique Fractures
- Spinal Fractures and Fusions
- Other Fusions
- Fractures stabilized by cast immobilization or surgically fixed with metal implantation.

The device is prescribed for patients who have a fracture that has presented with non-union or delayed union, or if they have risk factors that increase their chances of slow or non-healing. Risk factors can include fracture/injury complexity, previous medical conditions or habitual practices that have been shown to inhibit or slow the process of consolidation at the fracture site. Some examples of these factors include:

- Tobacco use and smoking
- Osteoporosis
- Diabetes
- Obesity
- Advanced Age
- Heart Conditions
- Currently on certain blood and/or non-steroidal anti-inflammatory drugs or medications
- Site Infection

Visual Indicators

The green LED light blinks rapidly throughout treatment. The red LED light located on the Response unit tells the user when their treatment time has elapsed. When the three hour treatment has concluded, the red LED light will blink continuously and the buzzer will beep three times. When the device has completed a total of 180 treatments, the LED light will blink five times and the buzzer will beep five times and it will not restart. Please contact Fintek if you continue to require treatment for non-healing at this time.

Frequently Asked Questions

Q: How does the device work?

A: The Bio-Electric Response uses Pulsed Electromagnetic Fields (PEMF) through capacitively coupled current to promote the healing of bone fractures. The device triggers a bone forming protein, serving as a cell receptor trigger to initiate the creation of new bone cells called osteoblasts. The osteoblasts are the building blocks of new bone and effectively fill in gaps that patients may have due to their fractures. As an additional benefit, PEMF may improve blood flow and vascularization to the treated area which potentially aids in healing tissues, ligaments, tendons, and nerves, while also reducing the risk of infection at the injury site.

Q: How do I use the Bio-Electric Response?

A: The device is connected to a set of leads and two electrodes. Place the electrodes on either side of your fracture site to create an arc. Turn the device on and treatment will begin. Fintek staff members are committed to ensuring that each patient is informed on how to properly use the device. In order to accomplish this, we ensure that we deliver the device in person wherever possible to show them how to properly use the device.

Q: How long do I need to wear the device for?

A: The device needs to be worn for 3 hours per day, every day. This continues until you are told by your physician that you have healed.

Q: Can I use the device for multiple treatments per day?

A: Yes you can. If you have two separate fracture locations you can wear the device on each location for three hours each.

Q: How long will the device last for?

A: The device is designed for 180 treatments. If you have completed 180 treatments and have not healed yet, we will either replace the device free of charge or re-code it for additional treatments.

Q: How long will the electrodes last?

A: The electrodes life time will vary, but they generally last 3-4 weeks or more if well maintained. We will provide you with 4 pairs of electrodes. Contact us if you need more.

Q: Will it be covered by my insurance?

A: Yes, but it does depend on the insurer and the information that is provided to them in order to support the claim. In our past experience, this treatment technology has been covered by WSIB/Worker's Compensation, Auto insurers, and quality supplementary health insurers. Fintek's office will

apply for this coverage on your behalf by submitting your signed claim form and assignment of benefits with a detailed claim letter and product efficacy documents and follow up with insurance companies for approval.

Q: How do I know the Bio-Electric Response is working?

A: You may start to see results on your clinic imaging within weeks of the beginning of treatment.

Q: Will I feel any pain?

A: Very few patients feel any sensation at all, but in some sensitive patients a mild increase in pain may be experienced in the first few days of use, this is due to the additional blood flow at the fracture site and generally goes away quickly.

Q: Is there any reason why I should not wear the device or are there any side effects?

A: You should not wear the device directly over or too near to a heart pace maker as it may interfere with the signal, and you should not wear it in many cases if you are pregnant. There are generally no side effects associated with wearing the device; it only heals bone that does not heal well and quickly. Some patients with very sensitive skin should consult their physician if any skin irritation occurs.

Q: What if I have any further questions?

A: You can contact the Fintek Bio-Electric Customer Service line and we will be able to assist you at any time.

Q: Is there a certain time of day that I should wear it?

A: No, you can wear the device at a time that suits you best. Ensure that you are wearing the device during a time that you're able to keep it on for a full three hours whenever possible.

Q: Do I need to charge the device?

A: No, the device runs on 9V batteries. The batteries will need to be changed when the red battery low light appears on the device. Fintek provides all patients with additional batteries for replacement.

Q: Can I split up the treatment time?

A: Fintek recommends that you wear the device for a full 3 hours per day, however, if you are required to take off the device, split up the treatment in two 1.5 hour treatments.

Limited Warranty

The Fintek Bio-Electric Response warranty covers any manufacturer defects for up to one-hundred and eighty treatments or for one year after the first treatment is delivered (whichever instance occurs first).

In cases where a defective Bio-Electric Response unit applies to the warranty, Fintek will replace the defective unit in a timely manner. In order to receive warranty service, please contact Fintek's toll-free telephone number at 1-866-440-5692 or email info@fintek.ca.

The Bio-Electric Response warranty does not apply to cases where defects or damage are caused by fire, flood, acts of god, or accidents, or if the device has been tampered with, misused or abused. In addition, the warranty is only valid for the patient that is prescribed the device, and cannot be transferred or assigned to any other individuals. Accessories used in concurrence with the device including the leads, electrodes, and batteries are not covered under this warranty.

Fintek holds the sole right to modify or extend the warranty only by written documentation. Physicians and sales representatives are not authorized to modify or extend the warranty.

Any defective products must be returned to Fintek Bio-Electric INC. to receive a replacement device.

This device is licensed solely for single patient use by the patient it is originally provided to. It must not be passed on or resold to others.

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